

Anti-NEU3 Antibody Picoband® Biotin Conjugated

Catalog Number: A04213-1-Biotin

About NEU3

This gene product belongs to a family of glycohydrolytic enzymes which remove sialic acid residues from glycoproteins and glycolipids. It is localized in the plasma membrane, and its activity is specific for gangliosides. It may play a role in modulating the ganglioside content of the lipid bilayer.

Overview

Product Name	Anti-NEU3 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UQ49

Technical Details

Immunogen	E.coli-derived human NEU3 recombinant protein (Position: H36-D414). Human NEU3 shares 68.8% and 67.7% amino acid (aa) sequence identity with mouse and rat NEU3, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NEU3 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.