

## Anti-MADCAM1 Antibody Picoband® Biotin Conjugated

Catalog Number: A04227-3-Biotin

### About MADCAM1

MADCAM1 (Mucosal Vascular Addressin Cell Adhesion Molecule 1), also known as MACAM1, is a protein that in humans is encoded by the MADCAM1 gene. By PCR-based analysis of somatic cell hybrids, this gene is mapped to chromosome 19. The protein encoded by this gene is an endothelial cell adhesion molecule that interacts preferentially with the leukocyte beta7 integrin LPAM-1 (alpha4 / beta7), L-selectin, and VLA-4 (alpha4 / beta1) on myeloid cells to direct leukocytes into mucosal and inflamed tissues. It is a member of the immunoglobulin superfamily and is similar to ICAM-1 and VCAM-1.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-MADCAM1 Antibody Picoband® Biotin Conjugated  |
| Reactive Species     | Human, Rat   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.                             |
| Host                 | Rabbit   |
| Uniprot ID           | Q13477   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthetic peptide corresponding to a sequence in the middle region of human MADCAM1, which shares 46.2% amino acid (aa) sequence identity with mouse MADCAM1. |
| Cross Reactivity    | No cross-reactivity with other proteins.  |
| Isotype             | Rabbit IgG  |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | Biotin  |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements.  |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MADCAM1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.