

Anti-p53DINP1 TP53INP1 Antibody

Catalog Number: A04229


Overview

Product Name	Anti-p53DINP1 TP53INP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p53DINP1 TP53INP1 Antibody catalog # A04229. Tested in WB,ICC,IF,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96A56


Technical Details

Immunogen	Full length recombinant protein of human p53DINP1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-p53DINP1 TP53INP1 Antibody (A04229) Images



Western blot analysis of p53DINP1 on different cell lysates using anti-p53DINP1 at 1/500 dilution. Positive control □ Line 1: SH-SY-5Y Line 2: HepG2 Line 3: Mouse stomach Line 4: Human stomach



ICC staining p53DINP1 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

2 Publications Citing This Product

1. PubMed ID: 10.3892/etm.2019.8007, Downregulation of microRNA-3934-5p induces apoptosis and inhibits the proliferation of neuroblastoma cells by targeting TP53INP1

2. PubMed ID: 33603839, Liu L,Zhang J,Liu Y. MicroRNA-1323 serves as a biomarker in gestational diabetes mellitus and aggravates high glucose-induced inhibition of trophoblast cell viability by suppressing TP53INP1. Exp Ther Med. 2021 Mar;21(3):230. doi:10.3892/etm.2021.9661.Epu

Visit [bosterbio.com/anti-p53dinp1-tp53inp1-antibody-a04229-boster.html](https://www.bosterbio.com/anti-p53dinp1-tp53inp1-antibody-a04229-boster.html) to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p53DINP1 TP53INP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.