

## Anti-ANO6 Antibody (N-term)

Catalog Number: A04251

### About ANO6

May act as a calcium-activated chloride channel.


### Overview

Product Name	Anti-ANO6 Antibody (N-term)
Reactive Species	Human
Description	Boster Bio Anti-ANO6 Antibody (N-term) (Catalog # A04251). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q4KMQ2

### Technical Details

Immunogen	This ANO6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 68-96 amino acids from the N-terminal region of human ANO6.
Predicted Reactive Species	Pig
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

## Anti-ANO6 Antibody (N-term) (A04251) Images



ANO6 Antibody (N-term) western blot analysis in SW480 cell line lysates (35ug/lane). This demonstrates the ANO6 antibody detected the ANO6 protein (arrow).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ANO6 Antibody (N-term)

For Research Use Only. Not for use in diagnostic procedures.