

Anti-MDH1 Antibody

Catalog Number: A04262

Overview

Product Name	Anti-MDH1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MDH1 Antibody catalog # A04262. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P40925

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human MDH1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-MDH1 Antibody (A04262) Images




Figure 1. Western blotting validation for Anti-MDH1 Antibody A04262

Western blot (WB) analysis of MDH1 polyclonal antibody at 1:500 dilution

Lane1:Hela cell lysate

Lane2:HepG2 cell lysate

Lane3:PC12 cell lysate

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MDH1 Antibody