

## Anti-ZW10 Antibody Picoband® HRP Conjugated

Catalog Number: A04266-1-HRP

### About ZW10

Centromere/kinetochore protein zw10 homolog is a protein that in humans is encoded by the ZW10 gene. This gene encodes a protein that is one of many involved in mechanisms to ensure proper chromosome segregation during cell division. This protein is an essential component of the mitotic checkpoint, which prevents cells from prematurely exiting mitosis.

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-ZW10 Antibody Picoband® HRP Conjugated   |
| Reactive Species     | Human   |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.    |
| Host                 | Rabbit  |
| Uniprot ID           | O43264  |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human ZW10 recombinant protein (Position: K35-K779).                                |
| Cross Reactivity    | No cross-reactivity with other proteins.   |
| Isotype             | Rabbit IgG   |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | HRP  |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ZW10 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.