

Anti-ING2 Antibody Picoband® Biotin Conjugated

Catalog Number: A04290-3-Biotin

About ING2

Inhibitor of growth protein 2 is a protein that in humans is encoded by the ING2 gene. This gene is a member of the inhibitor of growth (ING) family. Members of the ING family associate with and modulate the activity of histone acetyltransferase (HAT) and histone deacetylase (HDAC) complexes and function in DNA repair and apoptosis. Alternative splicing results in multiple transcript variants.

Overview

Product Name	Anti-ING2 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H160

Technical Details

Immunogen	E.coli-derived human ING2 recombinant protein (Position: Q6-K181). Human ING2 shares 96.6% amino acid (aa) sequence identity with mouse ING2.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ING2 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.