

## Anti-ABI2 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A04302-Fluoro550

### About ABI2

Abl interactor 2 also known as Abelson interactor 2 (Abi-2) is a protein that in humans is encoded by the ABI2 gene. Enables several functions, including SH3 domain binding activity; identical protein binding activity; and ubiquitin protein ligase binding activity. Contributes to small GTPase binding activity. Involved in Rac protein signal transduction; positive regulation of cellular component organization; and zonula adherens assembly. Acts upstream of or within peptidyl-tyrosine phosphorylation. Located in several cellular components, including filopodium tip; lamellipodium; and nucleoplasm. Part of SCAR complex. Is active in adherens junction. Colocalizes with actin filament.

### Overview

Product Name	Anti-ABI2 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NYB9

### Technical Details

Immunogen	E.coli-derived human ABI2 recombinant protein (Position: H151-D400).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ABI2 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.