

## Anti-PERQ2 GIGYF2 Antibody

Catalog Number: A04303-1

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-PERQ2 GIGYF2 Antibody   |
| Reactive Species     | Human, Mouse   |
| Description          | Boster Bio Anti-PERQ2 GIGYF2 Antibody catalog # A04303-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.        |
| Application          | ELISA, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.  |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | Q6Y7W6   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Synthesized peptide derived from human protein at AA range: 1071-1120   |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | PERQ2 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Suggested Dilutions | WB 1:500-2000<br>ELISA 1:5000-20000   |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PERQ2 GIGYF2 Antibody

For Research Use Only. Not for use in diagnostic procedures.