

## Anti-RGS7 Antibody Picoband® APC Conjugated

Catalog Number: A04305-1-APC

### About RGS7

Regulator of G-protein signaling 7 is a protein that in humans is encoded by the RGS7 gene. Enables G-protein beta-subunit binding activity and GTPase activator activity. Involved in positive regulation of GTPase activity. Located in cytosol.

### Overview

Product Name	Anti-RGS7 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P49802

### Technical Details

Immunogen	E.coli-derived human RGS7 recombinant protein (Position: M1-N483).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RGS7 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.