

Anti-GABA A Receptor Alpha 5 GABRA5 Antibody

Catalog Number: A04383-1


Overview

Product Name	Anti-GABA A Receptor Alpha 5 GABRA5 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-GABA A Receptor Alpha 5 GABRA5 Antibody catalog # A04383-1. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31644

Technical Details

Immunogen	Recombinant fusion protein of human GABRA5(NP_000801.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1000-1:4000

Anti-GABA A Receptor Alpha 5 GABRA5 Antibody (A04383-1) Images



Western blot analysis of extracts of various cell lines, using GABRA5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GABA A Receptor Alpha 5 GABRA5 Antibody

For Research Use Only. Not for use in diagnostic procedures.