

Anti-TPPP Antibody

Catalog Number: A04389


Overview

Product Name	Anti-TPPP Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TPPP Antibody catalog # A04389. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O94811

Technical Details

Immunogen	Recombinant fusion protein of human TPPP/p25(NP_008961.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-TPPP Antibody (A04389) Images



Western blot analysis of TPPP/p25 expression in Jurkat (A), mouse brain (B), rat brain (C) whole cell lysates.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TPPP Antibody

For Research Use Only. Not for use in diagnostic procedures.