

Anti-TPPP Antibody Picoband® Biotin Conjugated

Catalog Number: A04389-2-Biotin

About TPPP

Tubulin polymerization-promoting protein is a protein that in humans is encoded by the TPPP gene. Enables several functions, including GTPase activity; magnesium ion binding activity; and protein homodimerization activity. Involved in several processes, including microtubule cytoskeleton organization; negative regulation of tubulin deacetylation; and positive regulation of protein polymerization. Located in several cellular components, including mitochondrion; mitotic spindle; and perinuclear region of cytoplasm. Colocalizes with microtubule and microtubule bundle.

Overview

Product Name	Anti-TPPP Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O94811

Technical Details

Immunogen	E.coli-derived human TPPP recombinant protein (Position: M1-K219).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TPPP Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.