

## Anti-CDKN2D Antibody

Catalog Number: A04390

### Overview


Product Name	Anti-CDKN2D Antibody
Reactive Species	Human
Description	Boster Bio Anti-CDKN2D Antibody catalog # A04390. Tested in ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P55273

### Technical Details

Immunogen	Recombinant fusion protein of human CDKN2D(NP_001791.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	ICC/IF: 1:50-1:100

## Anti-CDKN2D Antibody (A04390) Images

---



Immunofluorescence analysis of NIH-3T3 cells using CDKN2D Rabbit pAb at dilution of 1:100 . Blue: DAPI for nuclear staining.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CDKN2D Antibody

For Research Use Only. Not for use in diagnostic procedures.