

Anti-Frizzled-7 FZD7 Antibody

Catalog Number: A04391


Overview

Product Name	Anti-Frizzled-7 FZD7 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Frizzled-7 FZD7 Antibody catalog # A04391. Tested in ELISA, IF, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, Flow Cytometry, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75084


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human FZD7. AA range:45-94
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 FC 1:50-200 IF/ICC1:100-500 ELISA 1:5000-20000


Anti-Frizzled-7 FZD7 Antibody (A04391) Images




Western Blot analysis of HELA cells using Frizzled-7 Polyclonal Antibody diluted at 1:2000



Western blot analysis of lysates from HT-29 cells, using FZD7 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of A549. 1, primary Antibody (red) was diluted at 1:200 (4°C overnight). 2, Goat Anti Rabbit IgG (H&L) - Alexa Fluor 594 Secondary antibody was diluted at 1:1000 (room temperature, 50min). 3, Picture B: DAPI (blue) 10min.



Immunofluorescence analysis of MCF7 cells, using FZD7 Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Frizzled-7 FZD7 Antibody

For Research Use Only. Not for use in diagnostic procedures.