

## Anti-Frizzled-7/FZD7 Antibody

Catalog Number: A04391-1

### Overview

Product Name	Anti-Frizzled-7/FZD7 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Frizzled-7/FZD7 Antibody catalog # A04391-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O75084

### Technical Details

Immunogen	E.coli-derived human FZD7 recombinant protein (Position: A32-W256).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 ELISA, 1:100-1000

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Frizzled-7/FZD7 Antibody

For Research Use Only. Not for use in diagnostic procedures.