

Anti-Frizzled-7 (Y87) FZD7 Antibody

Catalog Number: A04391Y87


Overview

Product Name	Anti-Frizzled-7 (Y87) FZD7 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Frizzled-7 (Y87) FZD7 Antibody catalog # A04391Y87. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75084


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 61-110 of Human Frizzled-7.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-Frizzled-7 (Y87) FZD7 Antibody (A04391Y87) Images



Western blot (WB) analysis of Frizzled-7 (Y87) pAb at 1:500 dilution Lane1:HeLa whole cell lysate(30ug) Lane2:BV2 whole cell lysate(30ug) Lane3:PC12 whole cell lysate(30ug)



Immunoprecipitation of HeLa cell lysates using Frizzled-7 pAb (Sepharose Bead Conjugate)#BD0048 (lane 2 and lane 3) and Nonspecific IgG Control (Sepharose Bead Conjugate)#BD0048 (lane 4 and lane 5) .Lane 1 is 30% input. The western blot was probed using Frizzled-7 pAb.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Frizzled-7 (Y87) FZD7 Antibody

For Research Use Only. Not for use in diagnostic procedures.