

## Anti-Snf1lk/SIK1 Antibody Picoband® HRP Conjugated

Catalog Number: A04399-2-HRP

### About SIK1

This gene encodes a serine/threonine protein kinase that contains a ubiquitin-associated (UBA) domain. The encoded protein is a member of the adenosine monophosphate-activated kinase (AMPK) subfamily of kinases that play a role in conserved signal transduction pathways. A mutation in this gene is associated with early infantile epileptic encephalopathy 30.

### Overview

Product Name	Anti-Snf1lk/SIK1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P57059

### Technical Details

Immunogen	E.coli-derived human Snf1lk/SIK1 recombinant protein (Position: M1-Q783). Human Snf1lk/SIK1 shares 80.5% and 81.2% amino acid (aa) sequence identity with mouse and rat Snf1lk/SIK1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Snf1lk/SIK1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.