

## Anti-PMS1 protein homolog 1 PMS1 Antibody

Catalog Number: A04425

### Overview


Product Name	Anti-PMS1 protein homolog 1 PMS1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PMS1 protein homolog 1 PMS1 Antibody catalog # A04425. Tested in IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P54277

### Technical Details

Form	Liquid
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB: 1:500~1:1000 IHC: 1:50~1:200  For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a>

### Anti-PMS1 protein homolog 1 PMS1 Antibody (A04425) Images

Figure 1. Western blotting validation for Anti-PMS1 protein homolog 1 PMS1 Antibody A04425



Western blot (WB) analysis of PMS1 polyclonal antibody at 1:500 dilution  
Lane1:Hela cell lysate  
Lane2:Raw264.7 cell lysate  
Lane3:PC12 cell lysate  
Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PMS1 protein homolog 1 PMS1 Antibody