

Anti-Bob 1 POU2AF1 Antibody

Catalog Number: A04431-1


Overview

Product Name	Anti-Bob 1 POU2AF1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Bob 1 POU2AF1 Antibody catalog # A04431-1. Tested in ELISA, IHC applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16633


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 201-250
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 ELISA 1:10000-20000 IF 1:50-200

Anti-Bob 1 POU2AF1 Antibody (A04431-1) Images



Immunohistochemical analysis of paraffin-embedded human colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Amygdala. 1, Antibody was diluted at 1:200 (4° overnight). 2, High-pressure and temperature EDTA, pH8.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 30min).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Bob 1 POU2AF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.