

Anti-IL-19 Antibody

Catalog Number: A04434


Overview

Product Name	Anti-IL-19 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-IL-19 Antibody catalog # A04434. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UHD0


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 51-100
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-2000 IHC 1:500-200 ELISA 1:10000-20000 IF 1:50-200

Anti-IL-19 Antibody (A04434) Images



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Western blot analysis of 293T HeLa VEC KB mouse-kidney lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-19 Antibody

For Research Use Only. Not for use in diagnostic procedures.