

## Anti-Cytochrome b5 (G82) CYB5A Antibody

Catalog Number: A04464-1

### Overview


Product Name	Anti-Cytochrome b5 (G82) CYB5A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cytochrome b5 (G82) CYB5A Antibody catalog # A04464-1. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P00167

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 51-100 of Human Cytochrome b5.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-Cytochrome b5 (G82) CYB5A Antibody (A04464-1) Images

---



Immunohistochemistry (IHC) analyzes of Cytochrome b5 (G82) pAb in paraffin-embedded human heart tissue.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome b5 (G82) CYB5A Antibody

For Research Use Only. Not for use in diagnostic procedures.