

## Anti-CD11a/ITGAL Antibody Picoband® Cy3 Conjugated

Catalog Number: A04466-1-Cy3

### About ITGAL

CD11a, also known as ITGAL or LFA-1, is a human gene which functions in the immune system. It is mapped to 16p11.2. CD11a is involved in cellular adhesion and costimulatory signaling. It is the target of the drug efalizumab. CD11a encodes the integrin alpha L chain. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This I-domain containing alpha integrin combines with the beta 2 chain (ITGB2) to form the integrin lymphocyte function-associated antigen-1 (LFA-1), which is expressed on all leukocytes. CD11a plays a central role in leukocyte intercellular adhesion through interactions with its ligands, ICAMs 1-3 (intercellular adhesion molecules 1 through 3), and also functions in lymphocyte costimulatory signaling. It is one of the two components, along with CD18, which form lymphocyte function-associated antigen-1.

### Overview

Product Name	Anti-CD11a/ITGAL Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P20701

### Technical Details

Immunogen	E.coli-derived human CD11a recombinant protein (Position: F161-L349). Human CD11a shares 73.7% amino acid (aa) sequence identity with mouse CD11a.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm

	Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

1. PubMed ID: 10.1016/j.ejphar.2005.12.080, Protective effect of bicyclol on acute hepatic failure induced by lipopolysaccharide and d-galactosamine in mice
2. PubMed ID: 25691929, Mu L, Jing C, Guo Z. Iran J Basic Med Sci. 2014 Nov;17(11):874-8. Expressions Of Cd11A, Cd11B, And Cd11C Integrin Proteins In Rats With Myocardial Hypertrophy.

Visit [bosterbio.com/anti-cd11a-picoband-trade-antibody-a04466-1-boster.html](http://bosterbio.com/anti-cd11a-picoband-trade-antibody-a04466-1-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD11a/ITGAL Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.