

Anti-PDC-E2 Monoclonal Antibody

Catalog Number: A04469


Overview

Product Name	Anti-PDC-E2 Monoclonal Antibody
Reactive Species	Human, Rabbit
Description	Boster Bio Anti-PDC-E2 Monoclonal Antibody catalog # A04469. Tested in WB applications. This antibody reacts with Human, Rabbit.
Application	WB
Clonality	Monoclonal 3F8
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P10515

Technical Details

Immunogen	Purified recombinant human PDC-E2 protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:1000-1:2000

Anti-PDC-E2 Monoclonal Antibody (A04469) Images



Western Blot analysis using PDC-E2 Monoclonal Antibody against Jurkat, A549, U251, F9, LNCAP, HeLa cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDC-E2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.