

Anti-Phosphoglycerate mutase 1 PGAM1 Antibody

Catalog Number: A04470-1


Overview

Product Name	Anti-Phosphoglycerate mutase 1 PGAM1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phosphoglycerate mutase 1 PGAM1 Antibody catalog # A04470-1. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P18669


Technical Details

Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:100 FC: 1:50-1:100

Anti-Phosphoglycerate mutase 1 PGAM1 Antibody (A04470-1) Images



Western blot analysis of PGAM1 on different lysates using anti-PGAM1 antibody at 1/2,000 dilution. Positive control:
Lane 1: A431 Lane 2: A549 Lane 3: Rat brain Lane 4: Mouse brain



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-PGAM1 antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phosphoglycerate mutase 1 PGAM1 Antibody

For Research Use Only. Not for use in diagnostic procedures.