

Anti-WDR79 WRAP53 Monoclonal Antibody

Catalog Number: A04484


Overview

Product Name	Anti-WDR79 WRAP53 Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-WDR79 WRAP53 Monoclonal Antibody catalog # A04484. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal 3F8
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9BUR4

Technical Details

Immunogen	Purified recombinant human WDR79 protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:1000-1:2000

Anti-WDR79 WRAP53 Monoclonal Antibody (A04484) Images



Western Blot analysis using WDR79 Monoclonal Antibody against HeLa, U2OS cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-WDR79 WRAP53 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.