

Anti-MMP-11 Antibody

Catalog Number: A04490


Overview

Product Name	Anti-MMP-11 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MMP-11 Antibody catalog # A04490. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P24347


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:500-1:1,000 IHC: 1:100-1:500

Anti-MMP-11 Antibody (A04490) Images



Western blot analysis of MMP11 on mouse spleen lysates using anti-MMP11 antibody at 1/1,000 dilution.



ICC staining of MMP11 in SW480 cells (green). The nuclear counterstain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

1 Publications Citing This Product

1. PubMed ID: 27884164, Ziwen Zhang, Detao Tao, Ping Zhang, Xue Liu, Yuchao Zhang, Jie Cheng, Hua Yuan, Laikui Liu and Hongbing Jiang Journal of Experimental & Clinical Cancer Research 2016;35:181 DOI: 10.1186/s13046-016-0458-0 Hyaluronan synthase 2 expressed by cancer-...

Visit bosterbio.com/anti-mmp-11-antibody-a04490-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MMP-11 Antibody

For Research Use Only. Not for use in diagnostic procedures.