

Anti-SIK2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A04515-1-Fluoro647

About SIK2

Enables ATP binding activity; magnesium ion binding activity; and protein kinase activity. Involved in intracellular signal transduction and protein autophosphorylation. Is active in endoplasmic reticulum membrane.

Overview

Product Name	Anti-SIK2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H0K1

Technical Details

Immunogen	E.coli-derived human SIK2 recombinant protein (Position: Q278-H919). Human SIK2 shares 88.2% amino acid (aa) sequence identity with mouse SIK2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SIK2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.