

Anti-MATE-1/SLC47A1 Antibody Picoband® Biotin Conjugated

Catalog Number: A04529-2-Biotin

About SLC47A1

Multidrug and toxin extrusion protein 1 (MATE1), also known as solute carrier family 47 member 1, is a protein that in humans is encoded by the SLC47A1 gene. This gene is located within the Smith-Magenis syndrome region on chromosome 17. It encodes a protein of unknown function.

Overview

Product Name	Anti-MATE-1/SLC47A1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96FL8

Technical Details

Immunogen	E.coli-derived human MATE-1/SLC47A1 recombinant protein (Position: M1-D508).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-MATE-1/SLC47A1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.