

## Anti-RBCK1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A04562-1-Fluoro550

### About RBCK1

Enables several functions, including identical protein binding activity; transcription coactivator activity; and ubiquitin binding activity. Involved in several processes, including positive regulation of canonical NF-kappaB signal transduction; protein linear polyubiquitination; and regulation of DNA-binding transcription factor activity. Part of LUBAC complex. Implicated in glycogen storage disease.

### Overview

Product Name	Anti-RBCK1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BYM8

### Technical Details

Immunogen	E.coli-derived human RBCK1 recombinant protein (Position: M1-R446). Human RBCK1 shares 90.3% and 89.9% amino acid (aa) sequence identity with mouse and rat RBCK1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-RBCK1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.