

## Anti-RBCK1 Antibody Picoband® HRP Conjugated

Catalog Number: A04562-1-HRP

### About RBCK1

Enables several functions, including identical protein binding activity; transcription coactivator activity; and ubiquitin binding activity. Involved in several processes, including positive regulation of canonical NF-kappaB signal transduction; protein linear polyubiquitination; and regulation of DNA-binding transcription factor activity. Part of LUBAC complex. Implicated in glycogen storage disease.

### Overview

Product Name	Anti-RBCK1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9BYM8

### Technical Details

Immunogen	E.coli-derived human RBCK1 recombinant protein (Position: M1-R446). Human RBCK1 shares 90.3% and 89.9% amino acid (aa) sequence identity with mouse and rat RBCK1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RBCK1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.