

Anti-CLEC9A Antibody HRP Conjugated

Catalog Number: A04577-HRP

About CLEC9A

C-type lectin domain family 9 member A is a protein that in humans is encoded by the CLEC9A gene. CLEC9A is a group V C-type lectin-like receptor (CTLR) that functions as an activation receptor and is expressed on myeloid lineage cells. By genomic sequence analysis, this gene is mapped to chromosome 12p13.31, centromeric to CLEC12B (617573) and telomeric to CLEC1A (606782).

Overview

Product Name	Anti-CLEC9A Antibody HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6UXN8

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human CLEC9A, which shares 47.5% and 55% amino acid (aa) sequence identity with mouse and rat CLEC9A, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-CLEC9A Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.