

Anti-Caldesmon (S783) CALD1 Antibody

Catalog Number: A04578S783


Overview

Product Name	Anti-Caldesmon (S783) CALD1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Caldesmon (S783) CALD1 Antibody catalog # A04578S783. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q05682


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 750-800 of Human Caldesmon.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-Caldesmon (S783) CALD1 Antibody (A04578S783) Images



Western blot (WB) analysis of Caldesmon (S783) polyclonal antibody at 1:500 dilution Lane1:MCF-7 whole cell lysate(20ug) Lane2:HeLa whole cell lysate(40ug) Lane3:A549 whole cell lysate(40ug) Lane4:PC12 whole cell lysate(40ug) Lane5:MEF whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of Caldesmon (S783) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caldesmon (S783) CALD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.