

Anti-PCMT1 Antibody

Catalog Number: A04579


Overview

Product Name	Anti-PCMT1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PCMT1 Antibody catalog # A04579. Tested in WB,ICC/IF,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P22061


Technical Details

Immunogen	Recombinant fusion protein of human PCMT1(NP_001238978.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:100 IP: 1:50-1:200


Anti-PCMT1 Antibody (A04579) Images



Western blot analysis of extracts of various cell lines, using PCMT1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunofluorescence analysis of HeLa cells using PCMT1 antibody. Blue: DAPI for nuclear staining.



Immunoprecipitation analysis of 200ug extracts of HepG2 cells using 1ug PCMT1 antibody. Western blot was performed from the immunoprecipitate using PCMT1 antibody at a dilution of 1:1000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PCMT1 Antibody

For Research Use Only. Not for use in diagnostic procedures.