

## Anti-Trx2 TXN2 Antibody

Catalog Number: A04586


### Overview

Product Name	Anti-Trx2 TXN2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Trx2 TXN2 Antibody catalog # A04586. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99757

### Technical Details

Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Trx2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-Trx2 TXN2 Antibody (A04586) Images



Western blot analysis of Trx2 expression in HEK293T (A), HeLa (B), mouse liver (C), rat lung (D), rat liver (E) whole cell lysates.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Trx2 TXN2 Antibody

For Research Use Only. Not for use in diagnostic procedures.