

Anti-RAB1B Antibody Picoband® HRP Conjugated

Catalog Number: A04589-1-HRP

About RAB1B

Ras-related protein Rab-1B is a protein that in humans is encoded by the RAB1B gene. Members of the RAB protein family, such as RAB1B, are low molecular mass monomeric GTPases localized on the cytoplasmic surfaces of distinct membrane-bound organelles. RAB1B functions in the early secretory pathway and is essential for vesicle transport between the endoplasmic reticulum (ER) and Golgi.

Overview

Product Name	Anti-RAB1B Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H0U4

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human RAB1B, which shares 100% and 96.8% amino acid (aa) sequence identity with mouse and rat RAB1B, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB1B Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.