

## Anti-CD42a/GP9 Antibody Picoband® HRP Conjugated

Catalog Number: A04595-2-HRP

### About GP9

Involved in blood coagulation and megakaryocyte development. Predicted to be part of glycoprotein Ib-IX-V complex. Is expressed in several structures, including brain; foregut; hemolymphoid system; liver; and male reproductive gland or organ. Human ortholog(s) of this gene implicated in Bernard-Soulier syndrome. Orthologous to human GP9 (glycoprotein IX platelet).

### Overview

Product Name	Anti-CD42a/GP9 Antibody Picoband® HRP Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O88186

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of mouse CD42a/GP9. Mouse CD42a/GP9 shares 89.5% amino acid (aa) sequence identity with human CD42a/GP9.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD42a/GP9 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.