

## Anti-Neurotensin Receptor 1 (NTS1) NTSR1 Antibody

Catalog Number: A04610-1

### Overview

Product Name	Anti-Neurotensin Receptor 1 (NTS1) NTSR1 Antibody
Reactive Species	Rat
Description	Boster Bio Anti-Neurotensin Receptor 1 (NTS1) NTSR1 Antibody catalog # A04610-1. Tested in IF, IHC, ICC applications. This antibody reacts with Rat.
Application	IF, IHC, ICC
Clonality	Polyclonal
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Chicken
Uniprot ID	P30989

### Technical Details

Immunogen	TVGTHNGLEHSTFNMTIEPGRVQAYC
Isotype	Chicken IgY
Form	Liquid
Purification	Affinity Purified
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

## Anti-Neurotensin Receptor 1 (NTS1) NTSR1 Antibody (A04610-1) Images




Figure 1. Immunofluorescent staining data of NTSR1 using Anti-Neurotensin Receptor 1 (NTS1) NTSR1 Antibody (A04610-1).

Immunofluorescent detection of NTS1 in cell bodies of neurons in rat dorsal root ganglia (red fluorescence).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Neurotensin Receptor 1 (NTS1) NTSR1 Antibody