

Anti-Neurotensin Receptor 1/NTSR1 Antibody Picoband® Biotin Conjugated

Catalog Number: A04610-4-Biotin

About NTSR1

Neurotensin receptor type 1 is a protein that in humans is encoded by the NTSR1 gene. Neurotensin receptor 1 belongs to the large superfamily of G-protein coupled receptors. NTSR1 mediates the multiple functions of neurotensin, such as hypotension, hyperglycemia, hypothermia, antinociception, and regulation of intestinal motility and secretion.

Overview

Product Name	Anti-Neurotensin Receptor 1/NTSR1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P30989

Technical Details

Immunogen	E.coli-derived human Neurotensin Receptor 1/NTSR1 recombinant protein (Position: A15-D400).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Neurotensin Receptor 1/NTSR1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.