

## Anti-ND3/MT-ND3 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A04611-1-Fluoro550

### About MT-ND3

Enables NADH dehydrogenase (ubiquinone) activity. Involved in mitochondrial electron transport, NADH to ubiquinone. Located in mitochondrial inner membrane. Part of respiratory chain complex I. Implicated in Leber hereditary optic neuropathy; Leigh disease; and Parkinson's disease.

### Overview

Product Name	Anti-ND3/MT-ND3 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P03897

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human ND3/MT-ND3.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-ND3/MT-ND3 Antibody - Fluoro550**

**For Research Use Only. Not for use in diagnostic procedures.**