

Anti-ND3/MT-ND3 Antibody Picoband® Cy3 Conjugated

Catalog Number: A04611-2-Cy3

About MT-ND3

Predicted to enable NADH dehydrogenase (ubiquinone) activity. Predicted to be involved in cellular response to glucocorticoid stimulus; mitochondrial electron transport, NADH to ubiquinone; and response to oxidative stress. Located in mitochondrion. Part of respiratory chain complex I. Is active in mitochondrial inner membrane. Is expressed in several structures, including alimentary system; brain; genitourinary system; hemolymphoid system gland; and liver and biliary system. Human ortholog(s) of this gene implicated in Leber hereditary optic neuropathy; Leigh disease; and Parkinson's disease. Orthologous to human MT-ND3 (mitochondrially encoded NADH:ubiquinone oxidoreductase core subunit 3).

Overview

| | |
|----------------------|--|
| Product Name | Anti-ND3/MT-ND3 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P03899 |

Technical Details

| | |
|---------------------|--|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of mouse ND3/MT-ND3. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ND3/MT-ND3 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.