

Anti-Neuroigin-3 NLGN3 Antibody

Catalog Number: A04627


Overview

Product Name	Anti-Neuroigin-3 NLGN3 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Neuroigin-3 NLGN3 Antibody catalog # A04627. Tested in WB applications. This antibody reacts with Human,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NZ94

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human NLGN3
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Neurologin-3 NLGN3 Antibody (A04627) Images



Western blot (WB) analysis of NLGN3 polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate Lane2:sp2/0 whole cell lysate

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Neurologin-3 NLGN3 Antibody

For Research Use Only. Not for use in diagnostic procedures.