

Anti-TRP 7/TRPC7 Antibody Picoband® Biotin Conjugated

Catalog Number: A04658-2-Biotin

About TRPC7

Transient receptor potential cation channel, subfamily C, member 7, also known as TRPC7, is a human gene encoding a protein of the same name. Using a cosmid/BAC contig, this gene is mapped to chromosome 21q22.3. TRPC7 is the most abundantly expressed long TRPC in peripheral blood and in some blood cell lines. The protein encoded by this gene is a calcium-permeable cation channel that is regulated by free intracellular ADP-ribose. The encoded protein is activated by oxidative stress and confers susceptibility to cell death.

Overview

Product Name	Anti-TRP 7/TRPC7 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9HCX4

Technical Details

Immunogen	E.coli-derived human TRP 7/TRPC7 recombinant protein (Position: R489-I862).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRP 7/TRPC7 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.