

## Anti-TRP 7/TRPC7 Antibody Picoband® FITC Conjugated

Catalog Number: A04658-2-FITC

### About TRPC7

Transient receptor potential cation channel, subfamily C, member 7, also known as TRPC7, is a human gene encoding a protein of the same name. Using a cosmid/BAC contig, this gene is mapped to chromosome 21q22.3. TRPC7 is the most abundantly expressed long TRPC in peripheral blood and in some blood cell lines. The protein encoded by this gene is a calcium-permeable cation channel that is regulated by free intracellular ADP-ribose. The encoded protein is activated by oxidative stress and confers susceptibility to cell death.

### Overview

Product Name	Anti-TRP 7/TRPC7 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9HCX4

### Technical Details

Immunogen	E.coli-derived human TRP 7/TRPC7 recombinant protein (Position: R489-I862).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRP 7/TRPC7 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.