

## Anti-EPHB3 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A04659-Fluoro488

### About EPHB3

Ephrin Receptor EphB3, is also known as human embryo kinase2(HEK2) or Eph-like tyrosine kinase2(ETK2). HEK2, which is a member of the EPH/ELK family of tyrosine kinases, encodes a 998-amino acid polypeptide having a single putative transmembrane domain, a secretory signal sequence, and 2 fibronectin repeats. The EPHB3 gene is mapped to human chromosome 3q21-qter. HEK2 interacts with 2 ligands of EPH-related kinases(LERKs), namely, LERK2(EFNB1) and LERK5(EFNB2). Coincubation of HEK2- and LERK2-expressing cells induces cell-cell adhesion and aggregation. Additionally, coexpression of HEK2 and LERK2 results in reduced kinase activity of HEK2.

### Overview

Product Name	Anti-EPHB3 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P54753

### Technical Details

Immunogen	E.coli-derived human EPHB3 recombinant protein (Position: H275-Q554).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EPHB3 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.