

## Anti-CARP (A243) ANKRD1 Antibody

Catalog Number: A04671-1


### Overview

Product Name	Anti-CARP (A243) ANKRD1 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-CARP (A243) ANKRD1 Antibody catalog # A04671-1. Tested in WB,IHC applications. This antibody reacts with Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15327

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 213-273 of Human CARP.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-CARP (A243) ANKRD1 Antibody (A04671-1) Images



Western blot (WB) analysis of ANKRD1 pAb at 1:500 dilution  
Lane1: The Muscle tissue lysate of Rat (40ug) Lane2: The Muscle tissue lysate of Mouse (40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CARP (A243) ANKRD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.