

## Anti-Lymphoid-specific helicase HELLS Antibody

Catalog Number: A04702


### Overview

Product Name	Anti-Lymphoid-specific helicase HELLS Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Lymphoid-specific helicase HELLS Antibody catalog # A04702. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NRZ9


### Technical Details

Immunogen	Recombinant full length Human HELLS.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:2000-1:5000 IHC: 1:50-1:200


## Anti-Lymphoid-specific helicase HELLS Antibody (A04702) Images



WesternBlot (WB) analysis of HELLS polyclonal antibody



Western blot analysis of extracts from normal (control) and HELLS knockout (KO) HeLa cells, using HELLS antibody.



Immunohistochemistry of paraffin-embedded mouse testis using HELLS Antibody.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Lymphoid-specific helicase HELLS Antibody

For Research Use Only. Not for use in diagnostic procedures.