

## Anti-TBX20 Antibody

Catalog Number: A04704

### Overview

Product Name	Anti-TBX20 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TBX20 Antibody catalog # A04704. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UMR3

### Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 211-260
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	TBX20 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

### 1 Publications Citing This Product

1. PubMed ID: , Single-Cell Profiles of Retinal Ganglion Cells Differing in Resilience to Injury Reveal Neuroprotective Genes

Visit [bosterbio.com/anti-tbx20-antibody-a04704-boster.html](https://www.bosterbio.com/anti-tbx20-antibody-a04704-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TBX20 Antibody

For Research Use Only. Not for use in diagnostic procedures.