

## Anti-RND3 Antibody

Catalog Number: A04727

### Overview

Product Name	Anti-RND3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RND3 Antibody catalog # A04727. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P61587

### Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 150-230
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	RND3 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

### 1 Publications Citing This Product

1. PubMed ID: -, Dong,H.,Sun,Q.,Zhang,Y.,Li,Y.,Yuan,F.,Mao,S., & Liu,B.(1899). Genetic deletion of Rnd3 suppresses apoptosis through NF $\kappa$ B signaling in the brain. Oncology Reports, 0, 0-0. <https://doi.org/10.3892/or.2020.7884>

Visit [bosterbio.com/anti-rnd3-antibody-a04727-boster.html](https://www.bosterbio.com/anti-rnd3-antibody-a04727-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-RND3 Antibody

For Research Use Only. Not for use in diagnostic procedures.