

Anti-DENND1A Antibody

Catalog Number: A04741-1

About DENND1A

DENND1A (DENN/MADD domain containing 1A), also known as Connecdenn, FAM31A or KIAA1608, is a 1,009 amino acid peripheral membrane protein. Recent studies has found DENND1A is strongly associated with PCOS (Polycystic ovary syndrome). DENND1A has several isoforms. This antibody(25658-1-AP) can recognize all the isoforms at 112 kDa, 52-64 kDa.

Overview

Product Name	Anti-DENND1A Antibody
Reactive Species	Human
Description	Boster Bio Anti-DENND1A Antibody catalog # A04741-1. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8TEH3

Technical Details

Immunogen	E.coli-derived human DENND1A recombinant protein (Position: D453-S557).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DENND1A Antibody

For Research Use Only. Not for use in diagnostic procedures.